



LABORATORY EQUIPMENT



WATER CHILLERS
THERMOSTATIC BATHS
DIGESTION SYSTEMS
ROTARY EVAPORATORS
VACUUM PUMPS
AIR COMPRESSORS
HOTPLATES
HOTPLATE STIRRERS



LabTech

WWW.LABTECHSRL.COM



THERMOSTATIC BATHS

LabTech Thermostatic Baths: flexibility and convenience

The LabTech cooling and heating circulators are the best choice for cooling small instruments. The LabTech circulators can be used in routine laboratory applications either with sample dipping or for external temperature control applications.

Features

- Temperature range -40°C ~ 100°C
- Both cooling and heating functions
- Internal and external liquid circulation
- High performance
- Ease of use
- Safe and quiet operation
- Robust and reliable
- P.I.D. precise temperature control
- High and low temperature protection
- Low liquid level protection
- Large LCD display
- Temperature control accuracy $\pm 0.2^{\circ}\text{C}$
- Excellent price/performance ratio



Main Applications

-
- | | |
|----------------------|-------------------|
| ■ Rotary evaporation | ■ Electrophoresis |
| ■ Distillation | ■ Polarimetry |
| ■ Calorimetry | ■ Photometry |
| ■ Fermentation | ■ Viscometry |
| ■ Refractometry | ■ Rheometry |
-

LabTech

WWW.LABTECHSRL.COM

LabTech, your reliable laboratory solution provider.

LabTech is a leading manufacturing company with skill and passion devoted to provide advanced laboratory solutions to the analytical community. We are located in Sorisole, northern Italy, with facilities in USA and China. Incorporated with global branches, LabTech is organized into R&D, production, marketing & sales and customer support. Major LabTech products include organic/inorganic sample preparation systems, laboratory cooling/heating/temperature control equipments, analytical instruments and gas generators. With knowledgeable scientists and experienced engineers, our R&D team has developed and released a series of new lab equipments and instruments dedicated to the modern analytical laboratories. With all of these accomplishments and our endless research and development efforts, we are proud to be your professional partner.



LabTech Srl

Via Fatebenefratelli, 1/5 - 24010 Sorisole (BG) - Italy
Tel: +39 035 576614 - Fax: +39 035 4729414
www.labtechsrl.com - email: marketing@labtechsrl.com