

TECHNICAL SPECIFICATIONS

SPECIFICATIONS	ED16S	ED36S	ED54S	EHD36S
Power (W)	1000	1400	3000	3000
Max Temperature (°C)	210	210	240	450
Accuracy (°C)	± 0.2	± 0.2	± 1	± 1
Tubes N°	16 (100mL)	36 (50mL)	54 (50mL)	36 (50mL)
Heating Block Material	Aluminium alloy PTFE coated	Aluminium alloy PTFE coated	Graphite PTFE coated	Graphite PTFE coated
Dimension mm (LxWxH)	390x380x190	420x420x190	536x398x290	538x442x198
Weight (Kg)	22	23	42	35

ACCESSORIES	GC-16	GC-36-L	GC-36-S	SC-36	TF-16	TF-36
Material	Glass	Glass	Glass	PP	PTFE	PTFE
Max Temperature (°C)	450	450	450	130	210	210
Height (mm)	120	180	120	105	125	125
Diameter (mm)	42	30	30	30	42	30
Volume (mL)	100	50	50	50	100	50

RACKS	BKS-16	BKS-36	BKT-36S	BKT-54	BKP-36S	BKP-54S
Material	Stainless steel	Stainless steel	PTFE Coated	PTFE Coated	Polycarbonate	PTFE
Dimension mm (LxWxH)	330x270x80	340x150x100	440x150x130	410x140x190	440x145x140	410x140x190
Diameter of holes (mm)	44	32	32	32	32	32
Holes quantity	16	18	18	18	18	18