

## OPTIONS FOR A COMPLETE PACKAGE

### VACUUM PUMPS

Totally oil free, mechanical vacuum pumps with features of advanced design, low noise level, high efficiency and long lifespan.

The LabTech vacuum pumps are specially designed for laboratory operations to meet the highest expectations of precision, reliability and ease of use. They can be used in a multitude of applications, such as vacuum filtration, degassing, impregnation, drying chambers, rotary evaporations, distillation, gel dryers, gas analytical systems and roughing of wide-range turbo molecular pumps.



	VP30	VP18R PLUS	VP50	VP50 PLUS
Power (W)	80	30	100	100
Pumping Speed (L/min)	30	18	50	50
Vacuum Pressure (mbar/psi)	250/3,62	20/0,29	150/2,17	150/2,17
Vacuum Controller	-	Manual	-	Manual
Noise (db)	≤50	≤42	≤58	≤58
Dimensions (LxWxH mm)	294x164x240	294x164x240	350x180x240	350x180x240
Weight (Kg)	3	3	6	6
Ordering information	LV5010	LV5080P	LV5050	LV5060